

# IOLMATIC®

—

## THE FUTURE IS PRELOADED



### Complete Versatility

- preloadable with any IOL design (hydrophobic & hydrophilic)
- fully- or semi-preloaded
- mini-incision  $\leq 2.2\text{mm}$

### Engineered for Perfection

- safe IOL delivery
- quick & easy handling
- reliable results

### Patented Technology

- hydrophobic IOL - Click & Shoot
- hydrophilic IOL - Double-Click & Shoot
- saves time in training, preparation & handling




©2014 IOLUTION  
Dammtorstraße 21b  
20354 Hamburg, Germany

Tel. +49 40 209347 100  
Fax +49 40 209347 109  
Mail [info@iolution.com](mailto:info@iolution.com)

# IOLMATIC® Universal Preloaded IOL System

## Hydrophobic IOL Click & Shoot

Preloaded IOLMATIC® in sterile blister



The folding element is pushed into the cartridge

*CLICK*  
the system is loaded



*SHOOT*  
haptics fold onto the IOL  
ensuring safe and reliable  
delivery



## Hydrophilic IOL Double-Click & Shoot

BSS-filled tray containing  
preloaded lens holder




IOLMATIC® in sterile blister

Lens holder preloaded  
with IOL



IOLMATIC® is  
connected to the  
lens holder




*1<sup>ST</sup> CLICK*  
IOLMATIC® is now  
preloaded

IOLMATIC® can remain  
inserted in the tray



The folding element is  
pushed into the cartridge

*2<sup>ND</sup> CLICK*  
the system is loaded



*SHOOT*  
haptics fold onto the IOL  
ensuring safe and reliable  
delivery



©2014 IOLUTION  
Dammtorstraße 21b  
20354 Hamburg, Germany

Tel. +49 40 209347 100  
Fax +49 40 209347 109  
Mail [info@iolution.com](mailto:info@iolution.com)